BD PureHub™ Disinfection Cap
Sterile
306598, 306599

BD Switzerland Sàrl
 Terre Bonne Park – A4
 Route de Crassier 17
 1262 Eysins, Switzerland
bd.com

TDS number: V201-002 – Rev. 01
 2020-June

1. General Information

1.1 Intended use

BD PureHub™ Disinfecting Cap are intended to be used as a disinfecting device for swabbable needle-free luer connectors prior to access and to act as a physical barrier between line accesses.

1.2 General description

BD PureHub™ Disinfecting Cap is comprised of a threaded plastic cap containing a porous pad pre-saturated with 70% Isopropyl Alcohol and sealed with a multilayer foil film.


BD PureHub™ Disinfecting Cap will disinfect the needle-free luer connector one (1) minute after application and act as a physical barrier for up to seven (7) days, if not removed. The effectiveness of BD PureHub™ Disinfecting Cap was tested in vitro against the following microorganisms: *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Escherichia coli*, *Pseudomonas aeruginosa*, *Acinetobacter baumannii*, *Candida glabrata* and *Candida albicans*. BD PureHub™ Disinfecting Cap reduced the number of microorganisms by greater than 4-log.


BD Catalog Number	BD Product Description	Compatibility when connected to other medical devices
306598	BD PureHub™ Disinfecting Cap Single	Female Luer connector of needleless connector
306599	BD PureHub™ Disinfecting Cap Strip	Female Luer connector of needleless connector

This document is BECTON DICKINSON property and cannot be reproduced or disclosed to any third party other than listed distribution list without BECTON DICKINSON written authorization. The content of this document may be subject to change without notice.

Bulk Configuration:

Strip configuration:

Further features: N/A

1.3 Certification

BD Catalog Number	BD Legal Manufacturer and ISO 13485 Certification	CE Certificate Number And Notified Body Brief Name	BD Manufacturing Site (Country of Origin) and ISO 13485 Certification	EC Representative (if applicable)
306598 306599	<p>Address: Becton Dickinson and Company 1 Becton Drive Franklin Lakes New Jersey 07417 United States</p> <p>ISO 13485 Certificate No.: MD19.2305</p>	<p>CE certified with NSAI Certificate No.: 252.1127</p>	<p>BD Medical Surgical 2153 12th Avenue Columbus NE 68601 United States</p> <p>ISO 13485 Certificate No.:MD19.2143</p>	<p>Becton Dickinson France SAS 11 rue Aristide Bergès ZI des Iles - BP 4 38801 Le-Pont-de-Claix France</p>

This document is BECTON DICKINSON property and cannot be reproduced or disclosed to any third party other than listed distribution list without BECTON DICKINSON written authorization. The content of this document may be subject to change without notice.

1.4 Materials

Component	Material
Cap	Polyethylene
Sponge	Polyester
IPA	The sponge is impregnated with USP grade 70% Isopropyl Alcohol (other names: 2-Propanol or Propan-2-ol)
Foil Lid (multilayer)	Polyester, aluminum foil and thermal adhesive

1.5 Materials of concern

Materials of concern are chemicals or substances that have been identified as having the potential to cause long term effects on humans or the environment.

Material	Comment
Phthalates	<p>Based on our ongoing data collection efforts and/or information received from our suppliers as per 26 May 2020, BD has not identified any</p> <ul style="list-style-type: none"> 1,2-Benzendicarboxylic acid, dihexyl ester (branched & linear) (CAS#68515-50-4), 1,2-Benzendicarboxylic acid, di-C6-8-branched alkyl esters (CAS#71888-89-6), 1,2-Benzendicarboxylic acid, di-C7-11-branched and linear alkyl esters (CAS#68515-42-4), 1,2-Benzendicarboxylic acid, di-C6-10 alkyl esters (CAS#68515-51-5), 1,2-Benzendicarboxylic acid, mixed decyl, hexyl, and octyl diesters (CAS#68648-93-1), Benzyl butyl phthalate (BBP) (CAS# 85-68-7), Bis (2-ethylhexyl) phthalate (DEHP) (CAS# 117-81-7), Bis (2-methoxyethyl) phthalate (DMEP) (CAS# 117-82-8), Di-n-hexyl phthalate (DnHP) (CAS# 84-75-3), Dibutyl phthalate (DBP) (CAS# 84-74-2), Diisobutyl phthalate (DIBP) (CAS# 84-69-5), Diisopentylphthalate (DIPP) (CAS# 605-50-5), Dipentyl phthalate (DPP) (CAS# 131-18-0), N-pentyl-isopentylphthalate (CAS# 776297-69-9) or Dicyclohexyl phthalate (DCHP) (CAS# 84-61-7) <p>in the articles and packaging with the product numbers as referenced above, in an individual concentration above 0.1% weight by weight (w/w). Furthermore, these products are not designed to contain any phthalates.</p>
Latex	Based on our ongoing data collection efforts and/or information received from our suppliers as per 26 May 2020, the articles with the product numbers referenced above are not formulated with natural rubber latex.
Bisphenol A	<p>Based on our ongoing data collection efforts and/or information received from our suppliers as per 26 May 2020, BD has not identified any</p> <ul style="list-style-type: none"> 4,4'-isopropylidenediphenol (BPA) (CAS# 80-05-7) <p>in the articles and packaging with the product numbers as referenced above. It is not a building block of any of the raw materials utilized and is not intentionally added. BD has not done any testing to evaluate levels of this chemical in these products.</p>
Substances of animal origin BSE/TSE	The raw materials used in the manufacture of these devices do not contain any animal tissue but may contain very small amounts of chemicals that are derived from animal-

This document is BECTON DICKINSON property and cannot be reproduced or disclosed to any third party other than listed distribution list without BECTON DICKINSON written authorization. The content of this document may be subject to change without notice.

	<p>derived raw materials. These products are manufactured using polymer resins which may contain very small amounts of stearic acid and related substances derived from tallow derivatives. Our resin suppliers have confirmed that these tallow-derived materials have been produced with multiple cycles of conditions at least as rigorous as those specified in Annex C.5 of EN 22442-1 and Section 6 of EMA 410/01 Rev. 3. Therefore, these raw materials meet or exceed the requirements of EN ISO 22442-1 and EMA 410/01 Rev. 3. Based on this information, these products are considered not to present any risk with respect to TSE/BSE or other animal-borne diseases. Furthermore, as recognized by MEDDEV 2.4/1, such derivative chemicals produced in accordance with the aforementioned standards and guidelines are considered irrelevant when determining the classification of a medical device (per MDD 93/42/EEC and EU No 722/2012).</p>
Polyvinyl chloride (PVC)	<p>The medical devices and packaging referenced above have not been designed nor intentionally manufactured with any additives or raw materials that contain PVC. No PVC is intentionally added to these medical products.</p>

1.6 REACH information

Based on our ongoing data collection efforts and/or information received from our suppliers as per 26 May 2020, BD has not identified any chemicals in the articles and packaging with the product numbers as referenced above, in an individual concentration above 0.1% weight by weight (w/w), which have been listed as SVHC and included in the "Candidate List" published by the European Chemical Agency (ECHA) on 15 January 2019 according to Art. 59 (1,10) of the Regulation (EC) N° 1907/2006 (REACH).

1.7 Biocompatibility

BD Medical products comply with the requirements of the standard for toxicity, pyrogenicity and biocompatibility of medical devices, ISO 10993 series - Biological Evaluation of Medical Devices.

1.8 Sterilization method

Sterilization method: **Gamma irradiation** as per EN ISO 11137 – "Sterilization of health care products – Radiation"

Guaranteed sterile until opening.

1.9 Shelf life and storage conditions

The BD PureHub™ shelf life has been assessed by stability studies in order to verify the functionality, physico-chemical and microbial properties over time.

BD PureHub™ references 306598 and 306599 have a shelf life of 3 years.

Store in dry and controlled temperature between 2 - 25°C and not exposed to strong light.

This document is BECTON DICKINSON property and cannot be reproduced or disclosed to any third party other than listed distribution list without BECTON DICKINSON written authorization. The content of this document may be subject to change without notice.

1.10 Standard

As per extract from the Declaration of Conformity:

Harmonized Standards	
EN 556-1:2001	Sterilization of Medical Devices - Requirements for Medical Devices to Be Designated "Sterile" - Part 1: Requirements for Terminally Sterilized Medical Devices
EN ISO 15223-1:2016	Medical Devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General Requirements
EN 1041:2008	Information supplied by the manufacturer of medical devices
EN 62366:2008	Medical devices. Application of usability engineering to medical devices
ISO 10993-1:2009	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
ISO 10993-4:2002 /amd.1:2009	Biological evaluation of medical devices - Part 4: Selection of tests for interactions with blood
ISO 10993-5:2009	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
ISO 10993-11:2009	Biological evaluation of medical devices - Part 11: Tests for systemic toxicity
ISO 10993-17:2009	Biological Evaluation of Medical Devices Part 17: Establishment of allowable limits for leachable substances
ISO 10993-18:2009	Biological evaluation of medical devices - Part 18: Chemical characterisation of materials
EN ISO 11137-1:2015	Sterilization of health care products - Radiation - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices
EN ISO 11137-2:2015	Sterilization of health care products - Radiation - Part 2: Establishing the sterilization dose
EN ISO 11737-1:2006	Sterilization of medical devices - Microbial methods- Part 1: Determination of a population of microorganisms on products
EN ISO 11737-2:2009	Sterilization of medical devices - Microbiological methods - Part 2: Tests of sterility performed in the validation of a sterilization process
ISO 11607-1:2009	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems
EN ISO 14971:2012	Medical devices. Application of risk management to medical devices.
EN ISO 13485:2016	Medical devices - Quality management systems - Requirements for regulatory purposes
EN ISO 14155:2011	Clinical investigation of medical devices for human subjects - Good clinical practice
Non-Harmonized Standards	
ISO 10993-10:2010	Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization
ISO 14644-1:2015	Cleanrooms and associated controlled environments - Part 1: Classification of air cleanliness by particle concentration
ISO TS 13004:2013	Sterilization of health care products - Radiation - Substantiation of selected sterilization dose: Method V _{DmaxSD}
ASTM D4169:2016	Standard Practice for Performance Testing of Shipping Containers and Systems
ASTM D3078-02:2013	Standard Test Method for Determination of Leaks in Flexible Packaging by Bubble Emission
ASTM F1980:2011	Standard Guide for Accelerated Aging of Sterile Barrier Systems for Medical Devices
ASTM F2825:2015	Standard Practice for Climatic Stressing of Packaging Systems for Single Parcel Delivery

Note:

The above standards reflect the status at the time of drafting this document. More information or updates are available on request in the Declaration of Conformity.

This document is BECTON DICKINSON property and cannot be reproduced or disclosed to any third party other than listed distribution list without BECTON DICKINSON written authorization. The content of this document may be subject to change without notice.

1.11 Classification

BD PureHub™ is classified as a **Class IIa** medical device under rule 15 of Annex IX of the Council Directive 93/42/EEC, as amended.

1.12 GMDN code

According to ISO 15225 (Medical devices - Quality management - Medical device nomenclature data structure), BD PureHub™ is referenced as follows:

GMDN Code: 61141

GMDN Term: Vascular catheter disinfection cap

1.13 Manufacturing practices

The entire manufacturing and testing processes are following the Manufacturing Practices as specified below:

- Incoming raw materials are verified via material inspection and testing and our suppliers are approved via our vendor management system.
- In addition to the automatic on-line inspections, in-process inspections are performed in addition to final product testing to ensure compliance with approved specifications.
- The manufacturing and testing details of each batch of product are recorded on a batch record which is retained in accordance with our document control procedures.
- BD operates a system of internal and external audits to maintain compliance.
- BD confirms that it will continue to adhere to relevant international standards in designing and manufacturing its products.

1.14 Other information

- (Material) Safety Data Sheets are not required for this product.
- Certificate of Food Contact (*Commission Regulation EU 1183/2012 on "plastic materials and articles intended for contact with food" and Directive 2002/72/CE (as amended) "relating to plastic materials and articles intended to come into contact with foodstuffs"*) is not required as BD products are used for general purpose injection and aspiration of fluids from vials, ampoules and parts of the body below the surface of the skin.
- Good Manufacturing Practices as defined by the FDA Pharmaceutical is not applicable for Medical Devices.
- Instruction for Use is available for these catalog numbers.

This document is BECTON DICKINSON property and cannot be reproduced or disclosed to any third party other than listed distribution list without BECTON DICKINSON written authorization. The content of this document may be subject to change without notice.

2. Packaging

2.1 Packaging configuration

BD Catalog Number	BD Product Description	Primary Packaging (Qty)	Shelf Box (Qty)	Shipping Case (Qty)	IFU Insert N/A / Yes / No*
306598	BD PureHub™ Disinfecting Cap Single	1	300	3000	Yes
306599	BD PureHub™ Disinfecting Cap Strip	1	300	4500	Yes

*"No": IFU may be available but not as an insert.


2.2 Packaging material

Component	Material
Unit Pack	Composed of White print mat, Polyester, Aluminum foil and Rayopeel thermal adhesive. The Peel Film is heat sealed to the Cap Housing flange
Shelf Box	Corrugated carton
Shipping Case	Corrugated carton
IFU	50 lb. white paper IFU is delivered as insert in shelf carton

2.3 Examples of labeling

Labels: According to European Medical Device directive, labels are multilingual (English, Spanish, Portuguese, French, German, Italian, Dutch, Swedish, Danish, Finnish, Greek and Norwegian)

Primary Packaging Label (Top Web) extracted from document 10000325454 related to reference 306598:


This document is BECTON DICKINSON property and cannot be reproduced or disclosed to any third party other than listed distribution list without BECTON DICKINSON written authorization. The content of this document may be subject to change without notice.

Primary Packaging Label (Top Web) extracted from document 10000294886 related to reference 306599:


Shelf Box extracted from document 10000269712 related to reference 306598:

This document is BECTON DICKINSON property and cannot be reproduced or disclosed to any third party other than listed distribution list without BECTON DICKINSON written authorization. The content of this document may be subject to change without notice.

Shipping Case extracted from document 10000269714 related to reference 306598:

Case Label extracted from document 10000295078 related to reference 306598:

This document is BECTON DICKINSON property and cannot be reproduced or disclosed to any third party other than listed distribution list without BECTON DICKINSON written authorization. The content of this document may be subject to change without notice.

IFU extract (English language) from document 10000269705 related to reference 306598 and 306599:

BD PureHub™
Disinfecting Cap (70% Isopropyl Alcohol)

BD PureHub™ Disinfecting Cap (70% Isopropyl Alcohol)
Tapón desinfectante BD PureHub™ (Alcohol isopropílico al 70%)
Tappe protettore BD PureHub™ (70% alcool isopropilico)
Capuchon de désinfection BD PureHub™ (Alcool isopropylique à 70%)
BD PureHub™ Desinfektionskappe (70 % Isopropylalkohol)
Cappuccio disinfettante BD PureHub™ (70% di alcol isopropilico)

1a.

1b.

2.

ENGLISH

Intended Use:
BD PureHub™ Disinfecting Caps are intended to be used as a disinfecting device for avoidable needle-free luer connectors prior to access and to act as a physical barrier between line accesses.

Description:
BD PureHub™ Disinfecting Cap is comprised of a three-sided plastic cap containing a porous pad pre-saturated with 70% Isopropyl Alcohol, and sealed with a multilayer foil film.

BD PureHub™ Disinfecting Cap will disinfect the needle-free luer connector one (1) minute after application and act as a physical barrier for up to seven (7) days, if not removed.

The effectiveness of BD PureHub™ Disinfecting Cap was tested in vitro against the following microorganisms: *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Escherichia coli*, *Pseudomonas aeruginosa*, *Acinetobacter baumannii*, *Candida glabrata* and *Candida albicans*. BD PureHub™ Disinfecting Cap reduced the number of microorganisms by greater than 4-log.

This product is not made with natural rubber latex or plasticizer Diethylhexyl Phthalate (DEHP).

CAUTIONS:

- Do not use if product is open or damaged.
- Verify the expiration date on the product label. Do not use if product has expired.
- Do Not Re-sterilize.
- After removing the foil film, visually inspect the contents of the cap and do not use if it appears dry or contains particulate matter.
- Do not place cap on a sterile field or aseptic area.
- Do Not Reuse. For Single Use Only. Re-use may lead to infection or other device injury.
- For use only on needle-free luer connectors.
- Always follow hospital protocol for disinfection of needle-free luer connectors.
- Store at controlled temperature (2-25°C)
- The safety of 70% Isopropyl Alcohol caps has not yet been directly established in clinical use scenarios with neuretics.

WARNING: Potential choking hazard.

WARNING: Keep out of reach of children. In case of accidental ingestion, seek medical help or consult a poison control center immediately.

WARNING: 70% Isopropyl Alcohol is not considered sporicidal and may not prevent Central Line-Associated Blood Stream Infection arising from bacterial spores (e.g., *Bacillus spp.*, *Clostridia*).

Instructions:

1. Peel the foil film from the BD PureHub™ Disinfecting Cap.
2. Apply cap to needle-free luer connector by pushing and seating.
3. Ensure that the cap is securely attached to the connector.
4. The BD PureHub™ Disinfecting Cap must remain on the needle-free luer connector for a minimum of one (1) minute and may remain on for up to seven (7) days.
5. Always place a new disinfecting cap on the needle-free luer connector after completing a line access procedure.
6. Discard cap after removal.

How Supplied:
BD PureHub™ Disinfecting Caps are provided in shelf boxes in the following configurations:
 Reorder No. 306598 Single (300 units)
 Reorder No. 306599 10 Shelf Strips (300 units/30 strips of 10)

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