

Selectavac Yankauers

Product Data Sheet

Description	Sterile, single use, surgical suction cannula with switch on and off valve, available with bulb tip and straight tip.
Intended Use	Selectavac Yankauers are used as a conduit to remove bodily fluids and secretions, surgical tissue debris and irrigation fluids. They are used in surgical sites and body cavities to facilitate observation and/or to clear an obstruction. (GMDN 35917)
Duration of Use	Normally intended for continuous use of less than 60 minutes, but will be patient specific depending on the duration of the surgical procedure.
Materials	Yankauer = Philips K Resin Handle = PS (Polystyrene) Suction Control Slider = PE
Product Information	Switch system allows the operator to manually open and close suction; Easy connection to Pennine range of tubing; Sterile, double wrapped, individual peel pouches; Latex Free.
Quality System	ISO 13485
EU Classification	Class IIa
Sterilised by	Ethylene Oxide
Shelf Life	5 years
Storage & Safety	Toxicity: Manufactured from biocompatible materials Storage: Store in a clean, cool, dry place away from direct sunlight



Product Code	Description	CH/Tip	Length	Qty per box
AEYS3009	Selectavac Suction, Bulb Tip with On/Off Slide Switch	24ch	265mm	50
AEYS3010	Selectavac Suction, Straight Tip with On/Off Slide Switch	24ch	265mm	50



GCMEDICA ENTERPRISE LTD., (WUXI)
Loujin Industrial Park, Shuofang, Wuxi 214143, Jiangsu, China



Luxus Lebenswelt GmbH
Kochstr.1, 47877, Willich, Germany

© 2021 Pennine Healthcare or one of its subsidiaries. All rights reserved.

The information given in this sheet is not intended to be exhaustive.

The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.