



Kidney Dish

FN-KD020101G01

- ✓ SINGLE USE
- ✓ MOULDED FIBRE PULP
- ✓ MACERATABLE
- ✓ MADE FROM RECYCLED PRODUCT

- Disposable single use product to reduce risk of cross infection
- Moulded products made from 100% recycled cellulose fibres.



Tested to Comply with
PAS29:1999

Regulatory Compliance:

- Compliant with: Medical Device Regulation 2017/745 Class 1 (non-sterile)
- Product Specification Assessment PAS 29:1999
- Disposability: By Maceration
- ISO 13485

Description:

- Capacity: 700ml
- Material: 100% Paper Pulp
- Dimensions: 250 x 140 x 45mm
- Colour: Grey

Quantity 300 per Carton

Outer Case Dimensions 455 x 260 x 310mm

Fluid Retention Testing

- Tested in Accordance with PAS 29:1999
- Guaranteed to hold fluid for 4 Hours
- Resistance Temperature: 35 °C (+/- 3)



Due to the nature of the manufacturing process the products may vary in appearance. This will not affect the performance of the product.

Finess Healthcare Group Co.,Ltd
No.1503 A1 Block, No.203 Lantian Road, Ningbo 315000 China
Tel:+0086-574-27835988 Fax:+0086-574-27835986