



Bedpan

FN-BP02010101



- ✓ SINGLE USE
- ✓ MOULDED FIBRE PULP
- ✓ MACERATABLE
- ✓ MADE FROM RECYCLED PRODUCT

- Disposable single use product to reduce risk of cross infection
- Moulded products made from 100% recycled cellulose fibres.

Regulatory Compliance:

- Product Specification Assessment PAS 29:1999
- Disposability: By Maceration
- ISO 13485

Description:

- Capacity: 2000ml
- Material: 100% Paper Pulp
- Dimensions: 360 x 300 x 85mm
- Colour: Grey

Fluid Retention Testing

- Tested in Accordance with PAS 29:1999
- Guaranteed to hold fluid for 4 Hours
- Resistance Temperature: 35 °C (+/- 3)

Quantity 100 per Carton

Outer Case Dimensions 510 x 380 x 310mm



Due to the nature of the manufacturing process the products may vary in appearance. This will not affect the performance of the product.

Finess Healthcare Group Co.,Ltd
No.1503 A1 Block, No.203 Lantian Road, Ningbo 315000 China
Tel:+0086-574-27835988 Fax:+0086-574-27835986