




REF: Grip-Lok® Medium Wide Adhesive (3300MWA, 3300MWA-NSB)

Description:

Grip-Lok® Medium Wide Adhesive is a disposable catheter stabilization device that locks PICC catheters in place. The Grip-Lok® Medium Wide Adhesive device can be placed on skin with gloved hands.

- Soft and flexible fabric design.
- Provides securement for both horizontal and vertical lifting, resulting from accidental line pulls.
- Simple to apply, inspect and adjust lines, tubes and catheters.
- Hypoallergenic, breathable and latex-free to reduce the risk of allergic reactions and skin irritation.

TECHNICAL SPECIFICATIONS
Materials
Bandage Adhesive: Urethane film laminated to a spun-laced non-woven tape, with a medical grade acrylic adhesive.
Standard Adhesive: Acrylic
Hook Fastener: Low-profile polypropylene
Loop: Low-profile knit nylon
Liner: Silicone-coated Kraft paper
Silicone Adhesive: A clear polyester film coated with a permanent silicone adhesive
Silicone Liner: Silicone fluorinated release liner
Other Features
Sterile: Individual, Single Use Packaging (3300MWA)
Non-Sterile: Bulk Packaging (3300MWA-NSB)
Not manufactured with natural Latex rubber
Hypoallergenic



Order Information

PN: 3300MWA	GTIN: 01618125134747	Unit of Purchase: 100/Box	Box Dimensions: 7.12 in x 6.12 in x 3.5 in (18.1 cm x 15.6 cm x 8.9 cm)
--------------------	-----------------------------	----------------------------------	---

Box Weight: 13.6 oz (385.6 g)	Expiration Date: 5 years from Date of Manufacture
--------------------------------------	--

PN: 3300MWA-NSB	GTIN: 00618125135110	Unit of Purchase: 2000/Box	Box Dimensions: 10.75 in x 8.75 in x 10 in (27.3 cm x 22.2 cm x 25.4 cm)
------------------------	-----------------------------	-----------------------------------	--

Box Weight: 8 lb-8 oz (3.9 kg)	Expiration Date: 5 years from Date of Manufacture
---------------------------------------	--



TIDI Products, LLC
570 Enterprise Drive
Neenah, WI 54956 USA
Toll Free: 800.521.1314
Phone: +1 920.751.4300
Fax: 800.837.7770

